




TRUMPF Med Australia

Operating Theatre Table Patient Positioning Accessories



TRUMPF

Med Aust

allen
amatech™

Allen Amatech Operating Theatre Table Patient Positioning Accessories

Elite and Standard Stirrups

Our Elite and Standard Stirrups represent the gold standard in lithotomy positioning. They allow for easy adjustment of abduction and lithotomy while maintaining the sterile field. The boot design reduces pressure under the popliteal fossa and the superficial peroneal nerve.

The stirrup enables safe and easy positioning while providing enhanced surgical site access



Arm and Hand Table Selection

We offer a hand table for every need. All of the tables are radiolucent with outstanding durability and are backed by a two-year warranty.

The carbon fiber, Carbon Lights™ tables have no frame or legs on the rectangular and minor procedure table to interfere with radiographic images.

The Basic tables utilize a phenolic board with a perimeter frame to allow for radiographic imaging.



LPS Arm Support

Safely and easily positions patients in lateral, prone and seated orientations. Allows for effective positioning of arms to prevent potential patient injury. The malleable hand rest enables surgeon to place patient's hand in any desired orientation.



Uro Catcher

Provides excellent fluid containment for cysto or hysteroscopy procedures on any table. Flexible drape support frame provides unobstructed surgeon access during manipulation and springs back to its original state to facilitate optimal fluid drainage.

Mounts easily with integral rail clamps to side rails and maintains its funnel shape.





Lift-Assist Beach Chair

Make Positioning Easy. Use Lift-Assist One-Handed Positioning. The Lift-Assist Beach Chair is the latest and most advanced option for positioning during shoulder procedures.

It offers unobstructed posterior access, and repositioning the patient is quick and easy.

The lift-assist design allows for nearly effortless positioning as the pistons support most of the patient's weight. Integral Clamps!



Easy Armboard with Trigger

With the trigger design, simply pull the radiolucent release mechanism for full articulation and release it to lock.

A spring lever accessible from both sides securely locks the armboard to the rail and eliminates the uncertainty of gravity locks. Eliminate Broken Armboard Teeth!



Deluxe Arthroscopic Legholder System

Smooth, uni-directional locking mechanism enables the caregiver to set rigid fixation with one hand operation. The hinged upper brace accommodates all size patients, while the structurally enhanced design enables the staff to apply pressure both medially and laterally.

The ratchet-style securement strap easily adjusts for precise amount of compression. Holds legs up to 89cm in circumference.



Radiolucent Spine Surgery Frame

Used to perform simple spine surgeries, has a patient weight capacity of 227 kg and is designed to fit standard operating theatre tables. This Frame enables the patient's abdomen to hang free during procedures. It is also radiolucent and allows for unrestricted C-arm access. The removable crank is used to adjust the amount of flex and tension in the patient supports. Every frame comes with a storage cart and one case of Skin Care Cover disposables.



Allen Amatech Operating Theatre Table Patient Positioning Accessories

Operating Theatre Table Clamps

We offer a variety of clamps for use with OT table accessories. The clamps are categorized by feature, function and rail type.

Clip-On Socket

The premium clamp instantly locks into place anywhere on the rail.



Tri-Clamp

Adjust the accessory while keeping the clamp locked securely on the rail.



Clip-On Blade Clamp

The premium clamp - instantly locks into place anywhere on the rail.



Rail Clamp

Basic design allows for securement of blade products.



Universal Clamp

Can be mounted vertically or horizontally on the OR table rail.



Simple Clamp

Can be sterilized and secured over the surgical drape.



Our product range includes

- Services Pendants
- LED Operating Theatre & Exam Lights
- HD Camera & Monitor Systems
- Operating Theatre Table Systems
- Operating Theatre Table Accessories
- ICU & Ward Hospital Beds
- Services Beam Systems
- Modular Room Systems
- Room Lighting
- Docking Trolley For Patient Transfer

TRUMPF



Med Aust

Unit 7, 9 McCabe Place, McCabe Estate CHATSWOOD NSW 2067
PO Box 309, WILLOUGHBY NSW 2068

P: (02) 9406 8900 • **F:** (02) 9406 8999 • **E:** info@trumpf-med.com.au
www.trumpf-med.com.au